
STANDARD Q COVID-19 Ag

STANDARD™ Q COVID-19 Ag Test **SD BIOSENSOR**


PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST

PREPARATION AND TEST PROCEDURE

■ PREPARATION

- Carefully read instructions for using the STANDARD Q COVID-19 Ag Test.
- Check the expiry date at the back of the foil pouch. Do not use the kit, if expiry date has passed.
- Check the test device and the desiccant pack in the foil pouch.


■ COLLECTION OF SPECIMEN

[Nasopharyngeal swab]


- Insert a sterile swab into the nostril of the patient, swab over the surface of the posterior nasopharynx. Withdraw the sterile swab from the nasal cavity.

- Insert the swab into an extraction buffer tube. While squeezing the buffer tube, stir the swab more than 5 times.


- Remove the swab while squeezing the sides of the tube to extract the liquid from the swab.

- Press the nozzle cap tightly onto the tube.


[Specimens in transport media]

- Apply a micropipette, collect the 350µl of specimen from the collection cup or VTM. Mix the specimen with an extraction buffer.
- Press the nozzle cap tightly onto the tube.

■ ANALYSIS OF SPECIMEN

- Apply 3 drops of extracted specimen to the specimen well of the test device.
- Read the test result in 15-30 minutes.

CAUTION Do not read test results after 30 minutes. It may give false results.

INTERPRETATION OF TEST RESULT

	* "C" Control Line	"T" Test Line
Negative		
Positive		
Invalid		

- A colored band will appear in the top section of the result window to show that the test is working properly. This band is control line (C).
- A colored band will appear in the lower section of the result window. This band is test line of SARS-CoV-2 antigen (T).
- Even if the control line is faint, or the test line isn't uniform, the test should be considered to be performed properly and the test result should be interpreted as a positive result.

* **The presence of any line no matter how faint the result is considered positive.**
 * **Positive results should be considered in conjunction with the clinical history and other data available.**

EXPLANATION AND SUMMARY

■ Introduction

Coronavirus is a single stranded positive-sense RNA virus with an envelope of about 80 to 120 nm in diameter. Its genetic material is the largest of all RNA viruses and is an important pathogen of many domestic animals, pets, and human diseases. It can cause a variety of acute and chronic diseases. Common signs of a person infected with a coronavirus include respiratory symptoms, fever, cough, shortness of breath, and dyspnea. In more severe cases, infection can cause pneumonia, severe acute respiratory syndrome, kidney failure, and even death. The 2019 new coronavirus, or "SARS-CoV-2 (COVID-19)", was discovered because of Wuhan Viral Pneumonia cases in 2019, and was named by the World Health Organization on January 12, 2020, confirming that it can cause colds and the Middle East Respiratory Syndrome (MERS) and more serious diseases such as acute respiratory syndrome (SARS). This kit is helpful for the auxiliary diagnosis of coronavirus infection. The test results are for clinical reference only and cannot be used as a basis for confirming or excluding cases alone.

■ Intended use

STANDARD Q COVID-19 Ag Test is a rapid chromatographic immunoassay for the qualitative detection of specific antigens to SARS-CoV-2 present in human nasopharynx. This test is for administration by healthcare workers and labs only, as an aid to early diagnosis of SARS-CoV-2 infection in patient with clinical symptoms with SARS-CoV-2 infection. It provides only an initial screening test result. This product is strictly for medical professional use only and not intended for personal use. The administration of the test and the interpretation of the results should be done by a trained health professional. The result of this test should not be the sole basis for the diagnosis; confirmatory testing is required.

■ Test principle

STANDARD Q COVID-19 Ag Test has two pre-coated lines, "C" Control line, "T" Test line on the surface of the nitrocellulose membrane. Both the control line and test line in the result window are not visible before applying any specimens. Mouse monoclonal anti-SARS-CoV-2 antibody is coated on the test line region and mouse monoclonal anti-Chicken IgY antibody is coated on the control line region. Mouse monoclonal anti-SARS-CoV-2 antibody conjugated with color particles are used as detectors for SARS-CoV-2 antigen device. During the test, SARS-CoV-2 antigen in the specimen interact with monoclonal anti-SARS-CoV-2 antibody conjugated with color particles making antigen-antibody color particle complex. This complex migrates on the membrane via capillary action until the test line, where it will be captured by the mouse monoclonal anti-SARS-CoV-2 antibody. A colored test line would be visible in the result window if SARS-CoV-2 antigens are present in the specimen. The intensity of colored test line will vary depending upon the amount of SARS-CoV-2 antigen present in the specimen. If SARS-CoV-2 antigens are not present in the specimen, then no color appears in the test line. The control line is used for procedural control, and should always appear if the test procedure is performed properly and the test reagents of the control line are working.

■ Kit contents

- Test device (individually in a foil pouch with desiccant)
- Extraction buffer tube
- Nozzle cap
- Sterile swab
- Instructions for use

KIT STORAGE AND STABILITY

Store the kit at 2-30°C / 36-86°F out of direct sunlight. Kit materials are stable until the expiration date printed on the outer box. Do not freeze the kit.

WARNINGS AND PRECAUTIONS

- Do not re-use the test kit.
- Do not use the test kit if the pouch is damaged or the seal is broken.
- Do not use the extraction buffer tube of another specimen.
- Do not smoke, drink or eat while handling specimen.

SPECIMEN COLLECTION AND PREPARATION

■ Transport medium

Virus Transport Medium(VTM)	Recommended Storage Condition	
	2°C to 8°C	25°C
Copan UTM™ Universal Transport Media	12 hours	8 hours
BD™ Universal Viral Transport	12 hours	8 hours
FA Transport Medium	12 hours	8 hours
ASAN PHARM UTM	12 hours	8 hours
GDL Korea UTM	12 hours	8 hours
Remel MicroTest M4RT	12 hours	8 hours
HBSS(+CaCl2)	12 hours	8 hours

NOTE When using viral transport medium (VTM), it is important to ensure that the VTM containing the specimen is warmed to room temperature. Cold specimens will not flow correctly and can lead to erroneous or invalid results. Several minutes will be required to bring a cold specimen to room temperature.

PERFORMANCE CHARACTERISTICS

■ Clinical evaluation

The sensitivity of the STANDARD Q COVID-19 Ag Test for rapid detection of SARS-CoV-2 antigen was established in prospective, multi institute, randomized, single-blinded studies conducted at a trial sites each in Brazil and India during the 2020 SARS-CoV-2 pandemic. A total of 115 positive specimens were tested using the STANDARD Q COVID-19 Ag Test. These specimens consisted of nasopharyngeal swabs from symptomatic patients. The specificity of STANDARD Q COVID-19 Ag Test was tested using 311 negative specimens. The sensitivity and specificity of the STANDARD Q COVID-19 Ag Test was compared to a commercialized molecular assay.

• Test sensitivity & specificity

The STANDARD Q COVID-19 Ag Test showed 96.52% of sensitivity and 99.68% of specificity.

- Summary of the sensitivity and specificity of the STANDARD Q COVID-19 Ag Test compared to PCR.

STANDARD Q COVID-19 Ag Test	PCR			Total
	Positive	Negative	Total	
	Positive	111	1	
Negative	4	310	314	
Total	115	311	426	
Sensitivity	96.52% (111/115, 95% CI 91.33 – 99.04%)			
Specificity	99.68% (310/311, 95% 98.22 – 99.99%)			

ANALYTICAL PERFORMANCE

- Limit of Detection (LoD): The study used "SARS-CoV-2 (2019-nCoV) NCCP 43326/2020 / Korea" strain. The titer of cultured virus was confirmed by PCR. The cell is inactivated and spiked into Nasopharyngeal swab specimen. The LoD is 1.25 X 10^{2.2} TCID₅₀/ml.

2019-nCoV Strain Tested		NCCP 43326/2020 / Korea									
Stock 2019-nCoV Titer		1 X 10 ^{2.2} TCID ₅₀ /ml									
Dilution	1/10	1/100	1/200	1/400	1/800	1/1600	1/3200	1/6400	1/12800		
Concentration in Dilution tested (TCID50/ml)	1 X 10 ^{2.2}	1 X 10 ^{2.1}	5 X 10 ^{2.1}	2.5 X 10 ^{2.1}	1.25 X 10 ^{2.1}	6.12 X 10 ^{2.1}	3.06 X 10 ^{2.1}	1.53 X 10 ^{2.1}	7.1 X 10 ^{2.1}		
Call rate of 5 replicates	100% (5/5)	100% (5/5)	100% (5/5)	100% (5/5)	100% (5/5)	80% (4/5)	0% (0/5)	0% (0/5)	0% (0/5)		
Call rates of 20 replicates near cut-off	NA	NA	NA	NA	100% (20/20)	56% (14/20)	0% (0/20)	NA	NA		
Lowest Concentration with Uniform Positivity per Analyte	1.25 X 10 ^{2.2} TCID ₅₀ /ml										
Limit of etection (LoD) per Virus Strain	1.25 X 10 ^{2.2} TCID ₅₀ /ml										

- Cross-Reactivity: There was no cross-reaction with potential cross-reactive substances except SARS-coronavirus.

Virus/Bacteria/Parasite	Strain	Source/Specimen type	Concentration	Results
SARS-coronavirus	Urbani	BEI/inactivated virus	3.5ug/ml	POS
	MERS-coronavirus	Jeddah_1_2013	Bionote/recombinant protein	10 ug/ml
Adenovirus	Type 1	Korea Bank for Pathogenic Viruses / live	3 X 10 ⁷ TCID ₅₀ /ml	NEG
	Type 3	National Culture Collection for pathogens / live	1.5 X 10 ⁴ TCID ₅₀ /ml	NEG
	Type 5	Korea Bank for Pathogenic Viruses / live	4 X 10 ⁷ TCID ₅₀ /ml	NEG
	Type 7	Korea Bank for Pathogenic Viruses / live	1.5 X 10 ⁴ TCID ₅₀ /ml	NEG
	Type 8	Korea Bank for Pathogenic Viruses / live	4 X 10 ⁷ TCID ₅₀ /ml	NEG
	Type 11	Korea Bank for Pathogenic Viruses / live	4 X 10 ⁷ TCID ₅₀ /ml	NEG
	Type 18	Korea Bank for Pathogenic Viruses / live	4 X 10 ⁷ TCID ₅₀ /ml	NEG
	Type 23	Korea Bank for Pathogenic Viruses / live	4 X 10 ⁷ TCID ₅₀ /ml	NEG
	Type 55	Korea Bank for Pathogenic Viruses / live	4 X 10 ⁷ TCID ₅₀ /ml	NEG
	H1N1 Denver	ATCC/live virus	3 X 10 ⁷ TCID ₅₀ /ml	NEG
	H1N1 WS/33	ATCC/live virus	3 X 10 ⁷ TCID ₅₀ /ml	NEG
Influenza A	H1N1 Pdm-09	ATCC/live virus	3 X 10 ⁷ TCID ₅₀ /ml	NEG
	H1N1 New Caledonia	ATCC/live virus	3 X 10 ⁷ TCID ₅₀ /ml	NEG
Influenza B	H1N1 New jersey	ATCC/live virus	3 X 10 ⁷ TCID ₅₀ /ml	NEG
	Nevada/03/2011	ATCC/live virus	3 X 10 ⁷ TCID ₅₀ /ml	NEG
Respiratory syncytial virus	B/Lee/40	ATCC/live virus	2.5 X 10 ⁴ TCID ₅₀ /ml	NEG
	B/Taiwan/2/62	ATCC/live virus	3 X 10 ⁷ TCID ₅₀ /ml	NEG
Respiratory syncytial virus	Type A	ATCC/live virus	3 X 10 ⁷ TCID ₅₀ /ml	NEG
Respiratory syncytial virus	Type B	ATCC/live virus	3 X 10 ⁷ TCID ₅₀ /ml	NEG
Legionella pneumophila	Bloomington-2	ATCC/live virus	5 X 10 ⁴ cells/ml	NEG
	Los Angeles-1	ATCC/live virus	5 X 10 ⁴ cells/ml	NEG
Mycobacterium tuberculosis	82A3105	ATCC/live virus	5 X 10 ⁴ cells/ml	NEG
	K		5 X 10 ⁴ cells/ml	NEG
	Erdman		5 X 10 ⁴ cells/ml	NEG
	HN878	Yonsel Univ. / inactivated and filter	5 X 10 ⁴ cells/ml	NEG
	CDC1551		5 X 10 ⁴ cells/ml	NEG
	H37Rv		5 X 10 ⁴ cells/ml	NEG
	4752-98 [Maryland (D1)6B-17]	ATCC/live	5 X 10 ⁴ cells/ml	NEG
	178 [Poland 23F-16]	ATCC/live	5 X 10 ⁴ cells/ml	NEG
	262 [CIP 104340]	ATCC/live	5 X 10 ⁴ cells/ml	NEG
	Slovakia 14-10 [29055]	ATCC/live	5 X 10 ⁴ cells/ml	NEG
Streptococcus pneumonia	Typing strain T1 [NCBI 11841, SF 130]	ATCC/live	5 X 10 ⁴ cells/ml	NEG
Streptococcus pyogenes	Mutant 22	ATCC/live	5 X 10 ⁴ cells/ml	NEG
Mycoplasma pneumoniae	FH strain of Eaton Agent [NCTC 10119]	ATCC/live	5 X 10 ⁴ cells/ml	NEG
Pooled human nasal wash - to represent diverse microbial flora in the human respiratory tract	M129-B7	ATCC/live	5 X 10 ⁴ cells/ml	NEG
	NA	Bionote / Normal pooled human nasal wash from healthy employees	NA	NEG
Coronavirus	229E	Zeptomatrix/inactivated	1 X 10 ⁴ TCID ₅₀ /ml	NEG
	OC43	Zeptomatrix/inactivated	1 X 10 ⁴ TCID ₅₀ /ml	NEG
	NL63	Zeptomatrix/inactivated	1 X 10 ⁴ TCID ₅₀ /ml	NEG

STANDARD Q COVID-19 Ag

STANDARD™ Q COVID-19 Ag Test **SD BIOSENSOR**

PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST

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