



CHLORAMPHENICOL EIA (5091CAP1p)


General

The Chloramphenicol (CAP) EIA is a competitive enzyme immunoassay for the detection of CAP and its major metabolite (CAP-glucuronide) in a variety of matrices. This kit has a unique sensitivity, robustness and reproducibility, which has resulted in a leading market position for already several years.

Kit characteristics

- Antibody cross reactivity:**

CAP	100%
CAP glucuronide	65%
CAP-base	< 1%
Thiamphenicol	< 1%
Florphenicol	< 1%



Chemical structure of chloramphenicol

- Kit contents:**

- Microtitre plate, 12 x 8, break 4 wells
- Sample dilution buffer (20 ml, 4x concentrated)
- Rinsing buffer (30 ml, 20x concentrated)
- Substrate solution (12 ml, Ready to use)
- Stop solution (15 ml, Ready to use)
- Conjugate (lyophilised)
- Antibody (lyophilised)
- CAP standards, 1 ml ready to use 0.025, 0.05, 0.1, 0.2, 0.5, 2 and 100 ng/ml
- Reconstitution/zero standard buffer (10 ml)


Assay procedure

Antibody, conjugate and sample/standard are pipetted into the wells of the microtitre plate and incubated for 1 hour at 2°C - 8°C. After a washing procedure ready to use substrate is added and incubated for 30 minutes at room temperature. Stop the reaction and read in a spectrophotometer at 450 nm. A calculation program is available upon request.

Inter Assay variation:	8%
Intra Assay variation:	6%

- Shelf life:**

Maximum 16 months.
Kit components are reusable after opening the kit



Assay characteristics

Matrices and sample preparation

- Urine, milk, milkpowder: direct after dilution, defatting.
- Urine, milk, serum: ethylacetate extraction, defatting.
- Tissue*¹, egg, honey: ethylacetate extraction, defatting.
- Feed: ethylacetate extraction, defatting.

LOD (ppb) ^{*2}	Recovery Rate ^{*2}
0.1	>90%
0.01	>90%
0.01/0.02	>80%
0.1	>90%

*¹ Tissues: meat, crab, shrimp, fish, liver.

*² LOD(Limit of Detection) and Recovery data; Validation according SANCO/1085/2000. Other matrices available upon request (Royal Jelly)