



CHLORAMPHENICOL FAST ELISA (5091CAPF)


General

The Chloramphenicol Fast (CAP) ELISA is a competitive enzyme immunosorbant assay for the detection of CAP and its major metabolite (CAP-glucuronide) in a variety of matrices. This kit has a unique sensitivity, robustness and reproducibility, which has resulted in a leading market position for already several years.

Kit characteristics

- Antibody cross reactivity:**

CAP	100%
CAP glucuronide	100%
Thiamphenicol	< 1%
Florphenicol	< 1%
- Kit contents:**
 - Microtitre plate, 12 x 8, break 4 wells
 - Sample dilution buffer (20 ml, 4x concentrated)
 - Rinsing buffer (30 ml, 20x concentrated)
 - Substrate solution (12 ml, Ready to use)
 - Stop solution (15 ml, Ready to use)
 - Conjugate (lyophilised)
 - Antibody (lyophilised)
 - CAP standards, 1 ml ready to use 0.04, 0.1, 0.2, 0.4, 1, 4 and 100 ng/ml
 - Reconstitution/zero standard buffer (10 ml)



Chemical structure of chloramphenicol

Assay procedure


Antibody, conjugate and sample/standard are pipetted into the wells of the microtitre plate and incubated for 30 minutes at room temperature (20°C - 25°C).

After a washing procedure ready to use substrate is added and incubated for 15 minutes at room temperature. Stop the reaction and read in a spectrophotometer at 450 nm.

A calculation program is available upon request.

Inter Assay variation: 8%
 Intra Assay variation: 6%

- Shelf life:**
 - Maximum 16 months.
 - Kit components are reusable after opening the kit



Assay characteristics

Matrices and sample preparation

- Urine
- Milk
- Serum
- Urine
- Milk
- Egg
- Tissue
- Liver
- Honey
- Feed

Procedure

- direct 0.5 ng/ml
- direct 0.2 ng/ml
- direct 0.2 ng/ml
- ethyl acetate extraction 0.02 ng/ml
- ethyl acetate extraction 0.02 ng/ml
- ethyl acetate extraction 0.02 ng/g
- ethyl acetate extraction 0.02 ng/g
- ethyl acetate extraction 0.02 ng/g
- ethyl acetate extraction 0.02 ng/g
- ethyl acetate extraction 0.5 ng/g

LOD (ppb)

LOD(Limit of Detection) and Recovery data; Validation according SANCO/1085/2000. Other matrices available upon request (Royal Jelly)