




PRESENTATION :
COMPASS[®] *LISTERIA METHOD*
FOR ENUMERATION



COMPASS[®] *Listeria*

Enumeration

- The **COMPASS[®] *Listeria*** method permits the enumeration of *L. monocytogenes* in human food products and environmental samples in as little as **2 days**.
- The **COMPASS[®] *Listeria* Agar** medium is based on the detection of 2 enzymatic activities :
 - β -D-Glucosidase specific to *Listeria* spp.
 - Phospholipase C, specific to *Listeria monocytogenes*.

COMPASS[®] *Listeria*

Enumeration

- This method is validated by AFNOR certification as compared to the reference method **NF EN ISO 11290-2/A1 (Feb. 2005)** : “Food Microbiology – Horizontal method for the detection and enumeration of *Listeria monocytogenes* - Part 1 : Enumeration method” and according to the standard **NF EN ISO 16140**.



BKR 23/05 - 12/07
Méthodes alternatives d'analyse
pour l'agroalimentaire
www.afnor-validation.org

- **COMPASS *Listeria*[®] Agar** is also validated for the detection of *Listeria* spp. and *Listeria monocytogenes* in human food products and environmental samples .
- **COMPASS *Listeria*[®] Agar** can also serve as the isolation medium used in the standard NF EN ISO 11290-2/A1.

COMPASS[®] *Listeria* Enumeration

Homogenization
of the sample



BPW

Primary dilution which can be used for additional analyses as referred to in NF EN ISO 6887-1

0.1 mL
streaking

OR

1 ml pour plate



Listeria spp. : Blue colonies without a halo
L. monocytogenes : Blue colonies with an opaque halo



T0

Test +

**CONFIRM' *L.mono*
broth**

SUSPENSION

INNOCULATION

48 ± 3 H à 37 ± 1 °C

CONFIRMATION

6 H à 37 ± 1 °C

COMPASS[®] *Listeria*

Enumeration

COMPASS[®] *Listeria* Agar

Excellent media performance & selectivity :

- Optimized formulation that favors *Listeria* growth.
- A judicious mixture of selective agents that insures inhibition of secondary flora.



COMPASS[®] *Listeria*

Enumeration

COMPASS[®] *Listeria* Agar

Ease of reading with clear distinctions :

- Blue colonies WITHOUT halos : *Listeria* spp.
- Blue colonies WITH opaque halos : *Listeria monocytogenes*.



COMPASS[®] *Listeria*

Enumeration

CONFIRM' *L.mono* broth

A rapid and visual confirmation test

- Easy reading : purple liquid turns yellow.
- Fast : results in 6 hours.

Reading the results can be made up to 24 hours later if needed.



Also exists in a solid media format : CONFIRM' *L.mono* Agar :

- Up to 6 confirmations per plate
- Results in 24 hours.

COMPASS[®] *Listeria*

Enumeration

The COMPASS[®] *Listeria* Agar enumeration method

Simple and rapid enumeration of *Listeria monocytogenes* :

- Negative results in only 2 days.
- No resuscitation step necessary.

COMPASS[®] *Listeria*

Enumeration

The COMPASS[®] *Listeria* Agar enumeration method

A flexible procedure :

- Enrichment and agar can be kept up to 3 days at 2-8 C before re-inoculations.
- Its formulation corresponds to that cited in the enumeration protocol for the standard EN ISO 11290-2 including the amendment A1.

SUSPENSION



INOCULATION

48 ± 3 H à 37 ± 1°C





CONFIRMATION

6 H à 37 ± 1 °C

COMPASS[®] *Listeria*

Enumeration



Listeria monocytogenes
and certain rare strains of *Listeria ivanovii*

Characteristic colonies :

Blue-green color surrounded by an
opaque halo.



***Listeria* spp.**

Characteristic colonies :

Blue-green color without any
halo.